



### **Product description**

Order number	97.8700.446
Product description	Charging station Multi Touch Pro, Charging station for Multi Touch Pro wireless scanner

### **Product characteristics**

Application	Obtaining blood and blood products
-------------	------------------------------------

### **Size**

Width of product	8 cm
Height of product	5 cm
Length of product	6 cm

### **Purity & certification**

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0039 cbm
Weight of product	0,35 kg
Weight of case	0,49 kg
EAN of inner box	4028339635720
EAN case	4028339635713

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

