



Product description

Order number	97.8700.449
Product description	Holder Multi Touch Pro, Table holder for Multi Touch Pro wireless scanner

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	10 cm
Height of product	6 cm
Length of product	30 cm

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0039 cbm
Weight of product	0.38 kg
Weight of case	0.52 kg
EAN of inner box	4028339635829
EAN case	4028339635812

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

