



Product description

| | |
|---------------------|--|
| Order number | 97.8700.461 |
| Product description | Wireless module DECT base station USB, for data communication between TOPSWING PRO II DECT or TRANSWAAG PRO and TRANSComm II PC software |

Product characteristics

| | |
|-------------|------------------------------------|
| Application | Obtaining blood and blood products |
|-------------|------------------------------------|

Size

| | |
|-------------------|--------|
| Width of product | 9 cm |
| Height of product | 4,5 cm |
| Length of product | 14 cm |

Purity & certification

| | |
|------------------|---|
| Product category | Appliance, according to machinery or equipment directive, CE |
| CE certified | CE |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

| | |
|----------------------|---------------|
| Minimum order qty. | 1 |
| Piece(s) / inner box | 1 |
| Piece(s) / case | 1 |
| Case volume | 0,0039 cbm |
| Weight of product | 0,16 kg |
| Weight of case | 0,3 kg |
| EAN of inner box | 4028339634327 |
| EAN case | 4028339634310 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

