



Product description

Order number	97.8700.461
Product description	Wireless module DECT base station USB, for data communication between TOPSWING PRO II DECT or TRANSWAAG PRO and TRANSComm II PC software

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	9 cm
Height of product	4.5 cm
Length of product	14 cm

Purity & certification

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0039 cbm
Weight of product	0.16 kg
Weight of case	0.3 kg
EAN of inner box	4028339634327
EAN case	4028339634310

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

