



Product description

Order number	97.8700.105
Product description	TRANSWAAG PRO, Blood mixing and weighing device, battery mode possible

Product characteristics

Application	Obtaining blood and blood products
Interface	1 x Anschluss für Scanner, 2 x USB, 1 x RS-232, 1 x RJ-45, 1 x Akku
Supply voltage	240 V
Power consumption during operation	24 W
Battery type	Lithium-ion battery
Battery power	6750 mAh

Size

Width of product	18 cm
Height of product	25 cm
Length of product	37.6 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product grey

Purity & certification

Product category Medical device
CE certified CE 0482 (MedCert)
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 42
Case volume 0.0341 cbm
Weight of product 2.01 kg
Weight of case 3.32 kg
EAN of inner box 4028339620122
EAN case 4028339620115

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

