



Product description

Order number	87.7708.909
Product description	Platelet agitator, for approx. 24 blood bags, 60 movements/minute (configurable), material: stainless steel, dimensions (LxWxH): 35 x 44.8 x 25 cm

Product characteristics

Application	Transfusion and infusion
Supply voltage	240 V
Power consumption during operation	50 W
Frequency	50 Hz - 60 Hz

Size

Width of product	44,8 cm
Height of product	25 cm
Length of product	35 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material stainless steel

Purity & certification

Product category no medical device | IVD

CE certified CE

Batched no

Packaging

Minimum order qty. 1

Piece(s) / inner box 1

Piece(s) / case 1

Piece(s) / pallet 8

Case volume 0,15 cbm

Weight of product 17 kg

Weight of case 19,5 kg

EAN of inner box 4038917140706

EAN case 4038917062169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

