



### **Product description**

Order number	87.7708.909
Product description	Platelet agitator, for approx. 24 blood bags, 60 movements/minute (configurable), material: stainless steel, dimensions (LxWxH): 35 x 44.8 x 25 cm

### **Product characteristics**

Application	Transfusion and infusion
Supply voltage	240 V
Power consumption during operation	50 W
Frequency	50 Hz - 60 Hz

### **Size**

Width of product	44,8 cm
Height of product	25 cm
Length of product	35 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material                      stainless steel

### **Purity & certification**

Product category                      no medical device | IVD

CE certified                              CE

Batched                                    no

### **Packaging**

Minimum order qty.                      1

Piece(s) / inner box                      1

Piece(s) / case                            1

Piece(s) / pallet                            8

Case volume                                0,15 cbm

Weight of product                        17 kg

Weight of case                            19,5 kg

EAN of inner box                        4038917140706

EAN case                                    4038917062169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

