



Product description

Order number	87.7708.916
Product description	Platelet incubator TI-2, 22 °C +/- 0.5 °C, CFC-free, material: stainless steel, door with double glazing, graphic data logger, monitoring of temperature, door position and agitator movement, dimensions (LxWxH): 64.8 x 104.4 x 140.3 cm

Product characteristics

Application	Transfusion and infusion
Supply voltage	240 V
Power consumption during operation	705 W
Frequency	50 Hz - 60 Hz

Size

Width of product	104.4 cm
Height of product	140.3 cm
Length of product	64.8 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material stainless steel

Purity & certification

Product category no medical device | IVD

CE certified CE

Batched no

Packaging

Minimum order qty. 1

Type of smallest
subpackaging case

Piece(s) / inner box 1

Piece (s) / outer case 1

Piece(s) / pallet 1

Case volume 1.35 cbm

Weight of product 156 kg

Weight of case 164.6 kg

EAN of inner box 4038917140775

EAN case 4038917062220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

