



Product description

Order number 87.7708.814
Product description Equipment rail for 4 MS 1000

Product characteristics

Application Obtaining blood and blood products

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 1
Case volume 0,0014 cbm
Weight of product 0,228 kg
Weight of case 0,42 kg
EAN of inner box 4038917341752
EAN case 4038917341745

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

