



Product description

Order number 87.7709.031
Product description Battery for AST-2 tube stripper

Product characteristics

Application Obtaining blood and blood products
Supply voltage 7.2 V
Battery power 700 mAh

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 432
Depth of box 41 mm
Width of box 41 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Height of box	102 mm
Depth of case	188 mm
Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.109 kg
Weight of case	0.26 kg
EAN of inner box	4038917141178
EAN case	4038917279697

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

