



Product description

Order number	87.7709.032
Product description	Battery charger for AST-2 tube stripper

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V
Power consumption during operation	17 W

Purity & certification

Product category	no medical device IVD
CE certified	CE supplier > no Sarstedt CE
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	432
Depth of case	188 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Width of case	188 mm
Height of case	100 mm
Case volume	0,0035 cbm
Weight of product	0,413 kg
Weight of case	0,56 kg
EAN of inner box	4038917141185
EAN case	4038917067157

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

