



Product description

Order number 87.7709.032
Product description Battery charger for AST-2 tube stripper

Product characteristics

Application Obtaining blood and blood products
Supply voltage 240 V
Power consumption during operation 17 W

Purity & certification

Product category no medical device | IVD
CE certified CE supplier > no Sarstedt CE
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 432
Depth of case 188 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification Charger

Page 2

| | |
|-------------------|---------------|
| Width of case | 188 mm |
| Height of case | 100 mm |
| Case volume | 0,0035 cbm |
| Weight of product | 0,413 kg |
| Weight of case | 0,56 kg |
| EAN of inner box | 4038917141185 |
| EAN case | 4038917067157 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT