



Product description

Order number 87.7709.033
Product description Mains adapter for AST-2 tube stripper

Product characteristics

Application Obtaining blood and blood products
Supply voltage 240 V

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD
CE certified CE supplier > no Sarstedt CE
Lot or serial number no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 432
Depth of case 188 mm
Width of case 188 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification
Mains adapter**

Page 2

Height of case	100 mm
Case volume	0,0035 cbm
Weight of product	0,365 kg
Weight of case	0,54 kg
EAN of inner box	4038917141192
EAN case	4038917067164

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-05-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT