



Product description

Order number	87.7709.035
Product description	Tube stripper, with auto-stop function, foot switch, table holder

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V
Power consumption during operation	20 W
Frequency	50 Hz - 60 Hz

Size

Width of product	140 mm
Height of product	71 mm
Length of product	250 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
CE certified	CE
Lot or serial number	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	84
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0,0174 cbm
Weight of product	2,103 kg
Weight of case	2,76 kg
EAN of inner box	4038917258890
EAN case	4038917258883

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

