



Product description

Order number	88.7708.095
Product description	Sealing clips, for sealing blood collection tubes, material: aluminium

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Material & colours

Product material	aluminium
Colour of product	grey

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1000
Piece(s) / inner box	1000
Piece(s) / case	99000
Piece(s) / pallet	6336000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0,0255 cbm
Weight of product	0,0001 kg
Weight of case	11,88 kg
EAN of inner box	4038917140454
EAN case	4038917550222

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

