



### **Product description**

Order number	97.8800.200
Product description	Stripping / Sealing Tongs, for stripping, segmenting and sealing blood collection tubes

### **Product characteristics**

Application	Obtaining blood and blood products
-------------	------------------------------------

### **Size**

Width of product	5 cm
Length of product	16 cm

### **Material & colours**

Colour of product	red
-------------------	-----

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## **Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0039 cbm
Weight of product	0,16 kg
Weight of case	0,34 kg
EAN of inner box	4028339660227
EAN case	4028339660210

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

