



Product description

Order number	97.8800.200
Product description	Stripping / Sealing Tongs, for stripping, segmenting and sealing blood collection tubes

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	5 cm
Length of product	16 cm

Material & colours

Colour of product	red
-------------------	-----

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0039 cbm
Weight of product	0,16 kg
Weight of case	0,34 kg
EAN of inner box	4028339660227
EAN case	4028339660210

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

