



### Product description

Order number	97.8750.156
Product description	ThermoScan USB adapter, for connecting the ThermoScan reading unit to the USB port of a PC

### Product characteristics

Application	Transport and storage
Interface	USB

### Size

Width of product	2 cm
Length of product	5.5 cm

### Purity & certification

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Packaging**

Minimum order qty.	1
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	1
Weight of product	0.0005 kg
EAN of inner box	4028339770421
EAN case	4028339770414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

