



## Product description

Order number	97.8750.000
Product description	TRANS Control 4 °C double element, for transporting erythrocyte concentrates and whole blood preparations at temperatures of 2 to 6 °C

## Product characteristics

Recommended temperature	4
Application	Transport and storage

## Size

Width of product	15.5 cm
Length of product	27 cm

## Material & colours

Colour of product	blue
-------------------	------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1  
Type of inner packaging                  case  
Piece(s) / inner box                      1  
Piece(s) / case                              1  
Case volume                                0.0039 cbm  
Weight of product                        0.85 kg  
Weight of case                              1.05 kg  
EAN of inner box                         4028339780123  
EAN case                                      4028339780116

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

