



Product description

Order number	97.8750.010
Product description	Temperature control element, double element, for transporting platelet concentrates at 22 °C

Product characteristics

Recommended temperature	22 °C
Application	Transport and storage

Size

Width of product	15,5 cm
Length of product	27 cm

Material & colours

Colour of product	green
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0,0039 cbm
Weight of product 0,85 kg
Weight of case 1,05 kg
EAN of inner box 4028339780321
EAN case 4028339780314

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

