



Product description

Order number	97.8750.016
Product description	Temperature control element, frame, for expanding the capacity of the TRANS Control 37°

Product characteristics

Recommended temperature	37 °C
Application	Transport and storage

Size

Width of product	15,5 cm
Length of product	27 cm

Material & colours

Colour of product	pink
-------------------	------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	yes

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0035 cbm
Weight of product	0,33 kg
Weight of case	0,53 kg
EAN of inner box	4028339780628
EAN case	4028339780611

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

