



Product description

Order number 97.8750.018
Product description TRANS Control -18 °C double element, for transporting frozen blood products

Product characteristics

Recommended temperature -18
Application Transport and storage

Size

Width of product 15.5 cm
Length of product 27 cm

Material & colours

Colour of product white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0.0035 cbm
Weight of product 1.2 kg
Weight of case 1.4 kg
EAN of inner box 4028339780727
EAN case 4028339780710

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

