



Product description

Order number	97.8750.020
Product description	TRANS Control -30 °C double element, for transporting FFP and frozen laboratory samples at temperatures from -30 to -32 °C

Product characteristics

Recommended temperature	-30
Application	Transport and storage

Size

Width of product	15.5 cm
Length of product	27 cm

Material & colours

Colour of product	yellow
-------------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 1
Case volume 0.0039 cbm
Weight of product 1.22 kg
Weight of case 1.42 kg
EAN of inner box 4028339780925
EAN case 4028339780918

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

