



Product description

| | |
|---------------------|--|
| Order number | 97.8750.021 |
| Product description | TRANS Control -30 °C frame, for expanding the capacity of the TRANS Control -30° |

Product characteristics

| | |
|-------------------------|-----------------------|
| Recommended temperature | -30 °C |
| Application | Transport and storage |

Size

| | |
|-------------------|---------|
| Width of product | 15,5 cm |
| Length of product | 27 cm |

Material & colours

| | |
|-------------------|--------|
| Colour of product | yellow |
|-------------------|--------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0,0039 cbm
Weight of product 0,29 kg
Weight of case 0,49 kg
EAN of inner box 4028339781021
EAN case 4028339781014

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

