



Product description

Order number	97.8750.142
Product description	Mono/Triple Box, for protecting multiple sensitive products against temperature changes, impact and shock, for 2 TRANS Control twin elements, with insulation, can be reduced to Mono Box or extended to Triple Box, material: styrofoam

Product characteristics

Application	Transport and storage
-------------	-----------------------

Size

Width of product	28.5 cm
Height of product	29 cm
Length of product	38.5 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material styrofoam
Colour of product grey

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 1
Weight of product 0.04 kg
Weight of case 0.74 kg
EAN of inner box 4028339781625
EAN case 4028339781618

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

