



Product description

Order number	87.7709.025
Product description	Tube stripper AST-2S, battery-powered, with auto-stop function, incl. transport case

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V
Power consumption during operation	20 W
Frequency	50 Hz - 60 Hz
Battery type	Nickel metal hydride battery

Size

Width of product	47 mm
Height of product	62 mm
Length of product	225 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	CE - manufacturer's self-declaration
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	84
Depth of box	388 mm
Width of box	358 mm
Height of box	125 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	1.896 kg
Weight of case	2.22 kg
EAN of inner box	4038917258869
EAN case	4038917258852

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

