



### **Product description**

Order number	87.7709.031
Product description	Battery Battery for AST-2 tube stripper

### **Product characteristics**

Application	Obtaining blood and blood products
Supply voltage	7.2 V
Battery power	700 mAh

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	yes

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Depth of box	41 mm
Width of box	41 mm
Height of box	102 mm
Depth of case	188 mm
Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.109 kg
Weight of case	0.24 kg
EAN of inner box	4038917141178
EAN case	4038917279697

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

