



Product description

Order number	87.7709.032
Product description	Charger for AST-2 tube stripper

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V
Power consumption during operation	17 W

Purity & certification

Product category	not regulated
CE certified	CE supplier
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification Charger

Page 2

Piece(s) / pallet	432
Depth of box	188 mm
Width of box	188 mm
Height of box	100 mm
Depth of case	188 mm
Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.413 kg
Weight of case	0.56 kg
EAN of inner box	4038917141185
EAN case	4038917067157

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT