



Product description

Order number	87.7709.033
Product description	Mains adapter Mains adapter for AST-2 tube stripper

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V

Purity & certification

Product category	no medical device IVD
CE certified	CE supplier
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	432
Depth of case	188 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Mains adapter

Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.365 kg
Weight of case	0.54 kg
EAN of inner box	4038917141192
EAN case	4038917067164

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

