



Product description

| | |
|---------------------|--|
| Order number | 87.7709.035 |
| Product description | Tube stripper AST-2F, with auto-stop function, foot switch, table holder |

Product characteristics

| | |
|------------------------------------|------------------------------------|
| Application | Obtaining blood and blood products |
| Supply voltage | 240 V |
| Power consumption during operation | 20 W |
| Frequency | 50 Hz - 60 Hz |

Size

| | |
|-------------------|--------|
| Width of product | 140 mm |
| Height of product | 71 mm |
| Length of product | 250 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| CE certified | CE |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 84 |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 125 mm |
| Case volume | 0.0174 cbm |
| Weight of product | 2.103 kg |
| Weight of case | 2.76 kg |
| EAN of inner box | 4038917258890 |
| EAN case | 4038917258883 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

