



### **Product description**

Order number	88.7708.095
Product description	Sealing clips, for sealing blood collection tubes, material: aluminium

### **Product characteristics**

Application	Obtaining blood and blood products
-------------	------------------------------------

### **Material & colours**

Product material	aluminium
Colour of product	grey

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	99000
Piece(s) / pallet	6336000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0001 kg
Weight of case	11.88 kg
EAN of inner box	4038917140454
EAN case	4038917550222

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

