



Product description

| | |
|---------------------|---|
| Order number | 74.5340.500 |
| Product description | TUR-BAG, Urine drainage system, 4 l, with pump ball, sterile, tube length: 100 cm, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------------|----------------|
| Colour of print/label | blue |
| Design | with pump ball |
| Application | Urine drainage |
| Volume of work | 4 l |

Size

| | |
|------------------|---------|
| Sample volume | 4 l |
| Max. volume | 4 l |
| Diameter of tube | 10.1 mm |
| Length of tube | 100 cm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-------------------------|
| Product material | Polyvinylchloride (PVC) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|---------------------|
| Product category | Medical device (MD) |
| CE certified | CE supplier |
| Purity standard | sterile |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 10 |
| Piece(s) / pallet | 560 |
| Depth of case | 395 mm |
| Width of case | 300 mm |
| Height of case | 235 mm |
| Case volume | 0.0278 cbm |
| Weight of product | 0.2455 kg |
| Weight of case | 2.46 kg |
| EAN of inner box | 6955824024395 |
| EAN case | 6955824024012 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

