



Product Specification
**TUR-BAG, Urine drainage system, 4 l, without
pump ball, sterile**

Page 2

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	960
Depth of case	395 mm
Width of case	295 mm
Height of case	260 mm
Case volume	0.0303 cbm
Weight of product	0.1645 kg
Weight of case	3.83 kg
EAN of inner box	6955824024401
EAN case	6955824024029

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP
support and is valid without signature.



SARSTEDT