




Product Specification
Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile

Page 1



Product description

Order number	74.5681.200
Product description	Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile, tube length: 120 cm, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Design	exchange bag 2,000 ml
Application	Urine drainage
Volume of work	500 ml
Anti-reflux valve	yes

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	120 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine measuring system, 2 l, with anti-reflux
valve, exchange bag 2,000 ml, sterile**

Page 2

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	10
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	280
Depth of box	400 mm
Width of box	300 mm
Height of box	55 mm
Depth of case	590 mm
Width of case	315 mm
Height of case	250 mm
Case volume	0.0465 cbm
Weight of product	0.3364 kg
Weight of case	4.04 kg
EAN of inner box	6955824024425
EAN case	6955824024043

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP
support and is valid without signature.



SARSTEDT