




**Product Specification**  
**Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 74.5681.200  |
| Product description | Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile, tube length: 120 cm, 1 piece(s)/blister |

**Product characteristics**

|                       |                       |
|-----------------------|-----------------------|
| Colour of print/label | blue                  |
| Design                | exchange bag 2,000 ml |
| Application           | Urine drainage        |
| Volume of work        | 500 ml                |
| Anti-reflux valve     | yes                   |

**Size**

|                  |         |
|------------------|---------|
| Sample volume    | 2 l     |
| Max. volume      | 2 l     |
| Diameter of tube | 10.1 mm |
| Length of tube   | 120 cm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine measuring system, 2 l, with anti-reflux  
valve, exchange bag 2,000 ml, sterile**

**Page 2**

**Material & colours**

|                   |                         |
|-------------------|-------------------------|
| Product material  | Polyvinylchloride (PVC) |
| Colour of product | transparent             |

**Purity & certification**

|                  |                   |
|------------------|-------------------|
| Product category | Medical device    |
| CE certified     | CE supplier       |
| Purity standard  | sterile           |
| Sterilisation    | Gas sterilisation |
| Batched          | yes               |

**Packaging**

|                                  |               |
|----------------------------------|---------------|
| Minimum order qty.               | 10            |
| Type of smallest<br>subpackaging | blister       |
| Piece(s) / inner box             | 1             |
| Piece (s) / outer case           | 10            |
| Piece(s) / pallet                | 280           |
| Depth of case                    | 580 mm        |
| Width of case                    | 290 mm        |
| Height of case                   | 235 mm        |
| Case volume                      | 0.0395 cbm    |
| Weight of product                | 0.3364 kg     |
| Weight of case                   | 4.04 kg       |
| EAN of inner box                 | 6955824024425 |
| EAN case                         | 6955824024043 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-27

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**