



Product description

| | |
|---------------------|--|
| Order number | 74.5684.002 |
| Product description | Urine measuring system exchange bag, 2 l, with drainage valve, sterile, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------------|---------------------|
| Colour of print/label | blue |
| Design | with drainage valve |
| Application | Urine drainage |
| Volume of work | 2 l |

Size

| | |
|---------------|-----|
| Sample volume | 2 l |
| Max. volume | 2 l |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-------------------------|
| Product material | Polyvinylchloride (PVC) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------|
| Product category | Medical device |
| CE certified | CE supplier |
| Purity standard | sterile |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|----------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 50 |
| Piece(s) / pallet | 4050 |
| Depth of case | 395 mm |
| Width of case | 245 mm |
| Height of case | 200 mm |
| Case volume | 0.0194 cbm |
| Weight of product | 0.0401 kg |
| Weight of case | 2.77 kg |
| EAN of inner box | 6955824024449 |
| EAN case | 6955824024067 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

