




**Product Specification**  
**Secretion bag, 1.5 l, tube length: 90 cm**

**Page 1**



**Product description**

Order number	74.5281.101
Product description	Secretion bag, 1.5 l, universal step connector, tube length: 90 cm, 200 piece(s)/case
PZN number	16333620

**Product characteristics**

Colour of print/label	blue
Design	universal step connector
Application	Drainage of exudate
Volume of work	1,5 l
Anti-reflux valve	no

**Size**

Sample volume	1.5 l
Max. volume	1.5 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

**Purity & certification**

Product category	Medical device
CE certified	CE supplier
Batched	yes

**Packaging**

Minimum order qty.	200
Type of smallest subpackaging	case
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	9600
Depth of case	580 mm
Width of case	390 mm
Height of case	150 mm
Case volume	0.0339 cbm
Weight of product	0.0325 kg
Weight of case	7 kg
EAN of inner box	6971243660945
EAN case	6971243660655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

