



Product description

Order number	97.8720.500
Product description	SAHARA 4 100 - 240 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 4 storage bags, operation via touchscreen, automatic aggregate status recognition by heating zone and ability to switch blood products conveniently while warming device is running

Product characteristics

Application	Warming
Power consumption during operation	1.000.00 W
Frequency	50 Hz - 60 Hz

Size

Width of product	574 mm
Height of product	348 mm
Length of product	554 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	Medical device equipment, CE
CE certified	CE - manufacturer's self-declaration
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of case	691 mm
Width of case	671 mm
Height of case	590 mm
Case volume	0.2736 cbm
Weight of product	23.7 kg
Weight of case	36.6 kg
EAN of inner box	4038917968737
EAN case	4038917968744

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

