



### **Product description**

Order number	14.1112
Product description	Membrane adapter, connection: Luer, for S-Monovette® needle/Multifly® needle, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free
PZN number	6796187

### **Product characteristics**

Example of use	Adaptation of an S-Monovette® Needle or (Safety) Multifly® Needle to a Luer system, e.g. Blood-Gas Monovette®
Suitable for	S-Monovette® needle/Multifly® needle
Connection	Luer

### **Size**

Length of product	16 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	translucent

### **Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of box	268 mm
Width of box	81 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	287 mm
Height of case	172 mm
Case volume	0.0184 cbm
Weight of product	0.0005 kg
Weight of case	1.49 kg
EAN of inner box	4038917085359
EAN case	4038917004053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

