



**Product Specification**  
**Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm**

**Page 1**



**Product description**

Order number	95.936
Product description	Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm, accessories for the Blood Gas Monovette® , 250 piece(s)/container

**Product characteristics**

Type of collection	capillary
Application	accessories for the Blood Gas Monovette®
Suitable for	blood gas capillaries

**Size**

Diameter	1 mm
Length of product	9 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm**

**Page 2**

**Material & colours**

Product material                      metal

**Purity & certification**

Product category                      In vitrodiagnostic device, CE

CE certified                              CE

Batched                                  yes

**Packaging**

Minimum order qty.                      250

Type of smallest  
subpackaging                              container

Piece(s) / inner box                      250

Piece (s) / outer case                      41250

Piece(s) / pallet                              412500

Depth of case                              268 mm

Width of case                                191 mm

Height of case                               100 mm

Case volume                                0.0051 cbm

Weight of product                           0.0001 kg

Weight of case                               2.43 kg

EAN of inner box                           4038917234610

EAN case                                    4038917234603

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**