



**Product Specification**  
**Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm**

**Page 1**



**Product description**

Order number	95.937
Product description	Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm, blood gas analysis, 10 piece(s)/case

**Product characteristics**

Type of collection	capillary
Application	blood gas analysis
Suitable for	blood gas capillaries

**Size**

Diameter	19 mm
Height of product	25 mm
Length of product	25 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product                      red

Purity & certification

Product category                      no medical device | IVD  
Batched                                    yes

Packaging

Minimum order qty.                    10  
Type of smallest subpackaging                      case  
Piece(s) / inner box                    10  
Piece (s) / outer case                   160  
Piece(s) / pallet                        51840  
Depth of case                            268 mm  
Width of case                            191 mm  
Height of case                            100 mm  
Case volume                              0.0051 cbm  
Weight of product                       0.0072 kg  
Weight of case                           1.29 kg  
EAN of inner box                        4038917234641  
EAN case                                  4038917234634

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

