



Product description

Order number	86.1996.062
Product description	ESR pipette, with printed graduation and assembled piston insert, with print

Product characteristics

Label/ Print	with print
--------------	------------

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of box	233 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR pipette, with print

Page 2

Width of box	164 mm
Height of box	76 mm
Depth of case	388 mm
Width of case	338 mm
Height of case	250 mm
Case volume	0.0328 cbm
Weight of product	0.0076 kg
Weight of case	8.61 kg
EAN of inner box	4038917135047
EAN case	4038917025966

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT