



Product description

Order number	90.1060.062
Product description	ESR rack, suitable for S-Monovette® , without back panel

Product characteristics

Design	without back panel
Suitable for	S-Monovette®

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1560
Depth of box	310 mm
Width of box	80 mm
Height of box	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
ESR rack, suitable for S-Monovette®

Page 2

Depth of case	430 mm
Width of case	325 mm
Height of case	130 mm
Case volume	0.0182 cbm
Weight of product	1.167 kg
Weight of case	23.8 kg
EAN of inner box	4038917124935
EAN case	4038917019545

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT