



**Product Specification**  
**ESR rack, suitable for S-Sedivette® , 10 measuring stations**

**Page 1**



**Product description**

Order number	90.1090
Product description	ESR rack, suitable for S-Sedivette® , 10 measuring stations

**Product characteristics**

Suitable for	S-Sedivette®
Number of positions	10

**Material & colours**

Product material	Polypropylene (PP)
------------------	--------------------

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	320
Depth of case	598 mm
Width of case	398 mm
Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.6018 kg
Weight of case	11.75 kg
EAN of inner box	4038917085304
EAN case	4038917004008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

