




**Product Specification**  
**ESR rack, suitable for Microvette<sup>®</sup> , 10 measuring stations**

**Page 1**



**Product description**

Order number 90.1091  
Product description ESR rack, suitable for Microvette<sup>®</sup> , 10 measuring stations

**Product characteristics**

Suitable for Microvette<sup>®</sup>  
Number of positions 10

**Material & colours**

Product material Polypropylene (PP)

**Purity & certification**

Product category In vitro diagnostic device, CE  
CE certified CE  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	8
Piece(s) / pallet	160
Depth of case	263 mm
Width of case	98 mm
Height of case	335 mm
Case volume	0.0086 cbm
Weight of product	0.2505 kg
Weight of case	5.1 kg
EAN of inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

