



Product description

| | |
|---------------------|---|
| Order number | 85.1019 |
| Product description | Safety lancet, Neonatal, pink, penetration depth: 1.2 mm, blade width: 1.5 mm, for capillary blood collection, for puncturing heel, finger, 200 piece(s)/case |

Product characteristics

| | |
|-------------------|--|
| Application | capillary blood collection |
| Example of use | Children, Adults when large blood volumes are required, neonates: standard blood count, blood sugar, bilirubin |
| Puncture site | heel, finger |
| Penetration depth | 1.2 mm |
| Safety product | yes |

Size

| | |
|-------------------|----------------|
| Blood volume | 50 µl - 250 µl |
| Width of product | 21 mm |
| Length of product | 52 mm |
| Width of blade | 1.5 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product pink

Purity & certification

| | |
|------------------|-------------------|
| Product category | Medical device |
| CE certified | CE 0124 (DEKRA) |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1200 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 200 |
| Piece (s) / outer case | 1200 |
| Piece(s) / pallet | 108000 |
| Depth of box | 118 mm |
| Width of box | 118 mm |
| Height of box | 151 mm |
| Depth of case | 366 mm |
| Width of case | 248 mm |
| Height of case | 173 mm |
| Case volume | 0.0157 cbm |
| Weight of product | 0.0029 kg |
| Weight of case | 4.08 kg |
| EAN of inner box | 4038917085069 |
| EAN case | 4038917057660 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

