



Product description

Order number	85.1026
Product description	Incision safety lancet, Safety-Heel® Newborn, light blue/white, penetration depth: 1 mm, blade width: 2.5 mm, for heel puncture, 50 piece(s)/case

Product characteristics

Example of use	Blood collection from newborn babies
Puncture site	Heel
Penetration depth	1 mm
Safety product	yes

Size

Width of product	35 mm
Length of product	40 mm
Width of blade	2.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Incision safety lancet, Safety-Heel® Newborn,
penetration depth: 1 mm

Page 2

Material & colours

Colour of product light blue/white

Purity & certification

Product category Medical device
CE certified CE supplier
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 200
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 200
Piece(s) / pallet 96000
Depth of case 273 mm
Width of case 262 mm
Height of case 160 mm
Case volume 0.0114 cbm
Weight of product 0.0094 kg
Weight of case 2.08 kg
EAN of inner box 6945630110474
EAN case 6945630110535

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT