



**Product description**

Order number	15.1670
Product description	Multivette® 600 Serum CAT, venous and capillary blood samples, preparation: Clotting activator, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap white, screw cap, flat base, 100 piece(s)/case

**Product characteristics**

Type of collection	Venous and capillary
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with print
Colour of print/label	white
Cap	standard screw cap
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base
Recommended centrifugation	10000 U/min
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Sample volume	600 µl
Diameter	8 mm
Length including cap	47.9 mm
Length excluding cap	44.4 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	201 mm
Width of box	81 mm
Height of box	81 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	178 mm
Case volume	0.016 cbm
Weight of product	0.0021 kg
Weight of case	2.76 kg
EAN of inner box	4038917087384
EAN case	4038917005814

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

