



### Product description

Order number	15.1673
Product description	Multivette® 600 Lithium heparin LH, venous and capillary blood samples, preparation: Lithium heparin, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap orange, screw cap, flat base, 100 piece(s)/case

### Product characteristics

Type of collection	Venous and capillary
Type of preparation	Lithium heparin
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Cap	standard screw cap
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base
Recommended centrifugation	2000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Multivette® 600 Lithium heparin LH, 600 µl, cap orange, screw cap**

**Page 2**

**Size**

Sample volume	600 µl
Diameter	8 mm
Length including cap	47.9 mm
Length excluding cap	44.4 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	201 mm
Width of box	81 mm
Height of box	81 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0022 kg
Weight of case	2.79 kg
EAN of inner box	4038917087414
EAN case	4038917005845

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**