



### Product description

Order number	15.1675.100
Product description	Multivette® 600 Lithium heparin gel LH, venous and capillary blood samples, preparation: Lithium heparin gel, nominal volume: 600 µl, (LxØ) with cap: 53.3 x 8 mm, with print, cap green, screw cap, flat base, 100 piece(s)/case

### Product characteristics

Type of collection	Venous and capillary
Type of preparation	Lithium heparin gel
Colour code	ISO
Label/ Print	with print
Colour of print/label	green
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base
Recommended centrifugation	10000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Sample volume	600 µl
Diameter	8 mm
Length including cap	53.3 mm
Length excluding cap	49.8 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	201 mm
Width of box	81 mm
Height of box	81 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0024 kg
Weight of case	3.01 kg
EAN of inner box	4038917133005
EAN case	4038917024648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

